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### Cation exchange capacity (CEC) – November 2013 – Louvain-la-Neuve

This is a measure of the negative charges as a function of the soil weight. This varies depending on the saturator cation, the salt concentration and the equilibrium pH.

The procedure can be summarized by the following steps:

- Soil saturation by a cation (it is adsorbed on the soil surface) and displacement of exchangeable cations
- Atomic absorption to measure the displaced cations concentrations
- Salt excess lavage with a solvent
- Saturator cation desorption with another salt
- Distillation of ammonium to titrate it and know its concentration


This method is applicable to all type of soils but in the case of calcareous soil a small amount of the carbonates is solved, rising up the value of the exchangeable Ca ( measured Ca = exchangeable Ca + portion or carbonated Ca).

#### Material required

- Balance and disposable plastic trays for weighting
- 5 g of brute soil sample, fraction < 2 mm
- Sand for percolation
- Ash less floc
- One column for each sample
- CEC trays
- One 250 mL volumetric flask for each sample
- Ammonium acetate
- 1 L Becker (it's easier to drop liquid on the probe from a Becker than from the 5L recipient)
- 5 L volumetric flask (to prepare ammonium acetate and KCl solutions), magnetic stirrer and stirring plate
- pH meter
- Ethanol 99%
- Plastic tray to collect ethanol
- Teflon
- 100 mL plastic flasks
- KCl
- Distillation setup, Erlenmeyer flasks, 100 mL distillation flasks, titration acids, magnetic stirrer, phenolphthalein, mixed indicator 4.

#### Extraction of exchangeable cations and preparation of ammonium solutions

##### 1° Preparation of the samples on columns

- 
- Rinse with water and Sand columns before using them
  - Put a small piece of Ash less the bottom of the column. Tight it a bit to ensure solid won't feel off but not too much, otherwise it will be difficult to filter the liquids. The function of this material is to be the support of the entire solid.
  - Put a small quantity of sand (a small spoon in our case) over the floc.
  - Weight some 5 g of soil and 20 g of sand. Mix it well and put over the first layer of sand.
  - Put a small quantity of sand on the top of column.
  - Prepare two samples without soil (20 g of sand) to measure if the sand is clean or if it has a supply of cations (Blanco).

2° Stick a piece of paper with the code of each sample written with pencil. Later, ethanol will be used and it could erase inscriptions made with a marker.

3° Place the columns on the CEC trays (they have two stages with holes to secure the position of the columns).

4° Prepare a 250 mL volumetric flask for each sample (clean and with the code written on it, in this case a marker can be used). Put each flask under the corresponding column.

5° Percolate 6 times 25 mL of ammonium acetate (1 N, pH 7) on each sample. Wait until all samples have percolated the liquid before starting the next round, this will make easier to keep the account and be sure that every sample had 6x25 mL of acetate. Collect the liquid on the flasks under the columns. Adjust the volume with milli-Q water, cover with a piece of Teflon and mix well. Analyse by atomic absorption in order to measure the exchangeable Ca, Mg, K and Na (place only 100 mL on a plastic flask, the rest is not necessary).

This step will displace the cations adsorbed on the soil and it will replace them with ammonium. It will be used an excess of ammonium to ensure the complete desorption of the soil cations and the complete adsorption of ammonium on all available positions.

***AVOID STOPPING FOR THE NIGHT WHILE THE ACETATE IS ON THE SOIL.***

6° Wash the excess of ammonium cations on the columns with ethanol 99%. It is not necessary to collect the liquid so a tray is put under the columns and the ethanol of all samples falls directly there. Empty on the sink.

*Fill completely the columns with ethanol and leave to percolate at least 5 times.*

This step is to wash the excess of ammonium that has not been adsorbed. All the adsorbed ammonium will remain in the soil.


7° Percolate 6 times 25 mL of KCl (10%, pH 3) following the same procedure than for the acetate. Similarly, collect the liquid on 250 mL volumetric flasks, adjust the volume and mix well. This liquid is being distilled to know the amount of ammonium that had been retained on the soil during the first step.

This step desorbs the ammonia retained on the soil, the  $NH_4^+$  are replaced by  $K^+$ . Knowing how much ammonia is recovered, it is possible to calculate the amount of negative charges of the soil.

#### *Dissolutions preparation*

- Ammonium acetate ( $CH_3COONH_4$ ) 1N: weight 77.08 g per litre (usually 5 L are prepared). If necessary, adjust the pH to 7 with a diluted solution of  $CH_3COOH$  or  $NH_4$ .
- Ethanol 99%: it is already manufactured like this.
- KCl 10% (in weight) pH 3: weight 100 g per litre (usually 5 L are prepared). If necessary, adjust the pH to 3 with a diluted solution of HCl or KOH.

#### *Ammonium distillation*

**1. Experimental set-up**


1. The system is heated with a gas flame. Initially the gas flow is set to the maximum to accelerate the heating process but once the water is already boiling it is lowered at the minimum needed to keep the boiling process.
2. The water flask is the vapour generator. The vapour is the mean to heat the sample and vaporize the ammonium. The water (demineralized) level has to be always between a minimum and a maximum, and it is always preferable to keep it near the maximum. The generated water vapour goes through the system and, if the vapour exit (4) is closed, it reaches the sample heating it.
3. This flask collects the vapour that is condensed in the system. This prevents this liquid to go back to the water flask, thing that would slow down the vapour generation process. It has to be emptied from time to time.
4. This point has to be closed when collecting the distillate. If this is open, the vapour will exit the system without reaching the sample. When changing the sample this has to be open so it will not be condensing vapour in the refrigerator.
5. This is the sample flask. Usually, 25 mL of sample are used but this depends on the volume of acid used during the titration (**if it is too small, inferior to 1 mL, it will be used more sample**). It has to be added a drop of phenolphthalein (to see the change to a basic medium) and one or two drops of sodium hydroxide (at the 50%) until the liquid changes to a pink colour. This is to ensure a basic medium that will convert all the ammonium into ammoniac (see reaction 1). It is placed on the set-up and secured with a spring. It is necessary that the liquid level reaches the extreme of the tube so vapour will enter the sample; to do so, demineralized water is added through the vapour exit (4). Then, the exit is closed.
6. The refrigerator will condense the vapour abandoning the sample. Water has to be circulating on the cooling jacket.
7. The liquid condensed on the refrigerator is collected on an Erlenmeyer flask. This flask has 10 mL of boric acid at the 5% in order to fix the distilled ammoniac (see reaction 2), and 2-3 drops of Mixed Indicator 5 in order to show the change of pH during the titration (initially, this liquid is purple and once the first drop of ammonium arrives, changing to a basic pH, in turns into green). Once the liquid starts to be collected, it is collected for 5 minutes. Only the first drops of liquid are ammonium, the rest is necessary to rinse the system and leave it clean for the next sample.

### *Solutions preparation*


H<sub>3</sub>BO<sub>3</sub>: 2% in weight. To prepare 1 L, it is necessary 20 g of solid boric acid). Heat while dissolving to facilitate the procedure.

### Titration and CEC calculation


Reactions involved on the distillation process:

*At the beginning of distillation, in the volumetric flask*

*When distillate starts to be collected, in the collector flask*

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This is to measure the quantity of ammonium that had been retained on the soil. While stirring with a magnetic bar, acid is added drop by drop until the colour changes from green to pink. This is the neutralization point, in this point it has been used the same amount of equivalents for the acid than for the neutralized base (in this case, H<sub>2</sub>BO<sub>3</sub><sup>-</sup>). Since reaction (2) indicates that the relation H<sub>2</sub>BO<sub>3</sub><sup>-</sup>/NH<sub>4</sub><sup>+</sup> is 1:1, the amount of equivalents for the ammonium will be the same than the amount for the H<sub>2</sub>BO<sub>3</sub><sup>-</sup> (they both form part of a salt) and this will allow the calculation of the equivalents of ammonium initially present in the 250 mL dissolution (the initial solution and the volume distilled have the same equivalents/L value). Knowing the equivalents/L in the 250 mL solution and the initial weight of soil, it will be possible to calculate the amount of ammonium retained on the soil for kg of soil.

At the neutralization point this equation can be applied:

$$(4) N_{acid} \cdot V_{acid} = N_{base} \cdot V_{base}$$
$$N_{\text{H}_2\text{BO}_3^-} = N_{\text{NH}_4^+} = \frac{N_{acid} \cdot V_{acid}}{V_{\text{NH}_4^+}} \left( \frac{\text{equiv}}{L} \right)$$

This value is the CEC and it is usually expressed in mili equivalents and referred to 100 gr of the dry matter. So, having the LOI calculated, we have to traduce the weight we took at the beginning of the procedure (brute weight) into dry matter and change the equivalents to mili equivalents and the gr to 100 gr.

$$N_{\text{NH}_4^+} \left( \frac{\frac{\text{equiv}}{L} \cdot 1000 \frac{\text{mequiv}}{\text{equiv}} \cdot 0.250 \text{ L (total volume of sample)}}{\text{g brute soil} \cdot \frac{\text{Dry matter \%} \left( \frac{\text{g dry soil}}{\text{g brute soil}} \right)}{100}} \cdot 100 \text{ g} \right) \equiv \text{mequiv} / 100 \text{ g}$$

Now, it is known the amount of mili equivalents that can be found in 100 g of soil.

**These units are equivalent to cmol-c/kg** (they take into account that not all cations contribute with the same number of charges). To pass from cmol/kg to cmol-c/kg it is necessary to just divide the cmol/kg (calculated with the molecular weight) between the number of charges of the cation (i.e. for Ca<sup>2+</sup> it would be divided by 2). If ph were 7, the CEC ( $\frac{\text{mequiv NH}_4^+}{100 \text{ g}}$ ) would be equal to the total addition of the exchangeable cations (Ca<sup>2+</sup>, K<sup>+</sup>, Mg<sup>2+</sup>, Na<sup>+</sup>).

The acid used depends on the quantity of ammonium present. Usually, we start with H<sub>2</sub>SO<sub>4</sub> 1/100 N because it is more precise than the other acid (HCl 1/20 N, more concentrated, thus, each drop has more protons and neutralizes faster the ammonium). If the necessary volume is too big it is necessary to switch to the more diluted acid. If the volume used is too small (inferior to 1 mL) more sample has to be used in the distillation.

It could be interesting to distillate the same sample 2 or 3 times in order to check the obtained value (this a subjective method and therefore the result could be wrong, by measuring it several times, accuracy is ensured).

*Use*

The CEC represents the negative positions available on the structure of clays to adsorb cations. The cations that are adsorbed there are called exchangeable cations because they can be easily replaced. The difference between the total addition of the exchangeable cations (S) and the CEC (T) indicates the saturation degree of the soil (V). It is easily calculated:

$$V = \frac{S}{T} \cdot 100$$

The CEC gets higher as the soil is weathered, the weathering of soil enhances clay formation, and clays are the minerals that have available negative positions where cations can be adsorbed. It is influenced also by the organic matter presence.

*Safety pictograms and remarks*

Ammonium acetate: use gloves and preferably also safety goggles

**Ethanol 99%:** use gloves and preferably also safety goggles. Be careful of vapours.

Signal word: Danger

Hazard statement(s)

- H225 Highly flammable liquid and vapour.

Precautionary statement(s)

- P210 Keep away from heat/sparks/open flames/hot surfaces. - No smoking.

**KCl:** wear gloves.

Hazard statement(s)

- H412 Harmful to aquatic life with long lasting effects.

Precautionary statement(s)

- P273 Avoid release to the environment.

**Boric acid:** use gloves and preferably also safety goggles

Signal word Danger

Hazard statement(s)

- H360FD May damage fertility. May damage the unborn child.

Precautionary statement(s)

- P201 Obtain special instructions before use.
- P308 + P313 IF exposed or concerned: Get medical advice/ attention.

**Sulphuric acid:** ALWAYS work on hood, use gloves and safety goggles. Rinse with plenty of water material used.

Signal word: Danger

Hazard statement(s)

- H314 Causes severe skin burns and eye damage.
- H412 Harmful to aquatic life with long lasting effects.

Precautionary statement(s)

- P273 Avoid release to the environment.
- P280 Wear protective gloves/ protective clothing/ eye protection/ face protection.
- P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
- P310 Immediately call a POISON CENTER or doctor/ physician.

**NaOH:** When making a solution with water pay attention to the heat released (exothermic reaction).

Signal word: Danger

Hazard statement(s)

- H290 May be corrosive to metals.
- H314 Causes severe skin burns and eye damage.

Precautionary statement(s)

- P280 Wear protective gloves/ protective clothing/ eye protection/ face protection.
- P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
- P310 Immediately call a POISON CENTER or doctor/ physician

**Phenolphthalein:** use gloves.

Signal word: Danger

Hazard statement(s)

- H341 Suspected of causing genetic defects.
- H350 May cause cancer.
- H361 Suspected of damaging fertility or the unborn child.

Precautionary statement(s)

- P201 Obtain special instructions before use.
- P281 Use personal protective equipment as required.
- P308 + P313 IF exposed or concerned: Get medical advice/ attention.

**Mixed indicator 5:**

Signal word: Danger

Hazard statement(s)

- H225 Highly flammable liquid and vapour.

Precautionary statement(s)


- P210 Keep away from heat/sparks/open flames/hot surfaces. - No smoking

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


Fig. 1 Relationship between soil to solution ratios and  $K_d$  at various percentages of adsorbed test substance

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
Fig. 2. Values of  $\log K_{ow}$  vs. pH for fluorescein (this study; Kasnavia et al., 1999) and eosin Y (this study). Data points represent average of 3 to 5 experiments, and error bars represent 1 standard deviation. Some error bars are smaller than the data symbol.

## 2 Octanol-water partition coefficients (Kow) vs. pH for fluorescent dye tracers (YASUHIRO OBA)

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
Soil type	pH range (in 0.01 M CaCl <sub>2</sub> )	Organic carbon content (%)	Clay content (%)	Soil texture*
1	4.5-5.5	1.0-2.0	65-80	clay
2	> 7.5	3.5-5.0	20-40	clay loam
3	5.5-7.0	1.5-3.0	15-25	silt loam
4	4.0-5.5	3.0-4.0	15-30	loam
5	< 4.0-6.0 <sup>s</sup>	< 0.5-1.5 <sup>s*</sup>	< 10-15 <sup>s</sup>	loamy sand
6	> 7.0	< 0.5-1.0 <sup>s*</sup>	40-65	clay loam/clay
7	< 4.5	> 10	< 10	sand/loamy sand

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